


COBRA-OS[®]

COBRA-OS[®]

Key Features



- 4 French profile
- Zone 1 marker (48 cm)
- Zone 3 marker (28 cm)
- Atraumatic flexible j-tip
- Balloon RO markers
- Suture tab
- 10cc Sword-handled syringe



INCLUDED _____ For Your Convenience



- 0.018" floppy tip guidewire
- 4Fr kink resistant sheath with excellent transition
- 21G 7cm beveled angiographic needle
- 10cc sword-handled syringe